

November 7, 2014

Welcome

Summary

- Cytoscape v3.2 is now available for immediate download.
- Major new features
 - Ability to show nodes as dynamic charts
 - RESTful interface enables “Cytoscape as a Service”
 - Faster startup, session load, and graph merge
 - Higher capacity for large graphs (both memory and draw speed)
 - Tentative support for Java 8
 - Export graph as a fully packaged web site ready for a web server
- Addresses 95 issues relative to v3.1.1 (see Section 10)
 - Various bugs and inconsistencies in the new Style pane
 - Minor issues in importing/exporting table data and PSIMI v2.5 XML files
 - Bugs in starting command scripts (via the Cytoscape command line)
 - Better user and App Developer documentation
 - Feature improvements for exporting networks to web pages
- 18 more Apps are available in the App Store as compared to the v3.1.1 release.
- Ready for use by all segments of the Cytoscape community.
- Continues v3.1.1’s core feature set, specifically:
 - Upward compatible and can coexist with previous Cytoscape versions
 - Most popular plugins from previous versions are available, with more to follow.
- What to do next ...

Release notes: <http://www.cytoscape.org/releasesnotes.html>

Dear Cytoscape User

We are proud and pleased to release Cytoscape v3.2. This new release contains both new features and improvements to functionality and documentation as described in Sections 7 and 10.

This release marks our most extensive testing efforts to date, including 51 different combinations of operating systems and Java Virtual Machines across Windows, Linux, and Mac OS platforms. Check out your system’s configuration here: <http://tinyurl.com/k5nppxo>. Note that v3.2 retains support for 32 bit JVMs and Windows XP, as they are a significant portion of Cytoscape installations. Note, too, that v3.2 introduces experimental support for Java 8.

We’re also pleased to report that in the last few months, Cytoscape downloads have averaged about 12,000 per month, with a peak of 15,000 in March. During each working day, Cytoscape is started over 2,000 times worldwide!

Finally, in v3.2, 18 Apps were added to the App Store (<http://apps.cytoscape.org>), bringing the total to 71. (So far in 2014, 39 apps have been added for v3.x, while only one plugin has been added for v2.x.) The apps new for v3.2 include:

ARACNE	CyToStruct	PanGIA	structureViz2
cddApp	EnrichmentMap	PetriScape	Synapse Client
cyLineUp	ENViz	PINBPA	WordCloud
CyNetworkBMA	iRegulon	Propagate	
cyREST	JGF App	RINalyzer	

We expect to make a v3.2.1 maintenance release in mid Q1 2015, followed by a v3.3 feature release a few months later. As our plans firm up, we'll publish them to the Cytoscape community. For more on our roadmap, see Section 9.

*If you are a new or casual Cytoscape user, **welcome aboard!** Feel free to either read on or jump ahead to *What to Do Next* for a quick start.*

For clarity, references to "Cytoscape 3.x" apply to all Cytoscape versions 3.0.0 and later.

For Cytoscape v2.x Users

Previous Cytoscape releases (through v2.8.3) focused on adding numerous annotation, import, layout, analysis, and visualization features, all of which have combined to deliver critical insights into biological processes of interest to Cytoscape's broad and loyal user community. In fact, many of these features were created and contributed by Cytoscape users, themselves, using Cytoscape's plugin architecture.

Version 3.0 marked a major reorganization of Cytoscape program code – it enabled us to deliver additional features more quickly and reliably than before, including improvements that benefit new users, casual users, power users, and app developers. Along the way, some features were changed to allow them to better interoperate with other features, some little-used features were dropped, and user interface was improved. We were particularly excited about the new Cytoscape App Store (described below), which afforded users access to the growing collection of *apps* (formerly called *plugins*) produced by fellow Cytoscape users. A table of new features is presented in Section 7.

Since v3.0, we released updates v3.0.1 and v3.0.2 as maintenance releases, and then v3.1 and v3.2 as feature releases, including greatly improved visual styles, node filtering, speed, network capacity, and web publishing.

While all version 3 releases are upwardly compatible with previous Cytoscape versions (i.e., they can load, analyze, and display networks and attributes stored using previous versions, including v2.x), they can coexist with previous versions on the same workstation. It is not necessary to uninstall previous versions before installing a version 3 release, and you can choose to use both on the same workstation.

Version 3 is rapidly maturing and has achieved a reputation for speed, reliability, and support that makes it a *superior choice* (compared to v2.x) for most uses. If you are using v2.x, now may be a good time to switch to v3.x and enjoy the new Apps available in the App Store (<http://apps.cytoscape.org>).

A table describing Cytoscape support is presented in Section 8.

We are very excited to make this release, and hope you find it an indispensable tool in your research. To get you started, we have provided additional introductory material to get you started on the right foot, and we have also provided background material for general interest.

Please let us know how we're doing and how we can better advance your research!

The Cytoscape Team

What to Do Next

We hope to help you use Cytoscape in the most productive way. Your next steps depend on how experienced you are with Cytoscape and what you would like to do. In the table below, please find the category that best describes you, and proceed accordingly.

Type of Cytoscape User	Where to Read Next
New	For New Users (Section 1)
Casual	For Casual Users (Section 2)
Power	For Power Users (Section 3)
App Developers	For App Developers (Section 4)
General Interest	Try getting oriented by reading: General Background (Section 5) Cytoscape App Store (Section 6) Then try browsing: New Features (Section 7) Support Offering (Section 8) Future Plans (Section 9) To get your feet wet, try: For New Users (Section 1)

1. For New Cytoscape Users

What is Cytoscape?

Cytoscape is software that can help you visualize complex networks (as *graphs*) representing relationships between genes, proteins, or other entities. Cytoscape tracks and displays these entities (as graph *nodes*) and relationships between them (as graph *edges*). Cytoscape, per se, does not assign meaning to either nodes or edges – you can assign your own meanings to them, and you can add your own node or edge attributes (e.g., name, measurements, states, etc.) that customize your graph to represent semantics important to you. While Cytoscape provides simple analytics and visualizations, others are available as apps from the Cytoscape App store (see Section 6).

What Do I Need Before I Can Use Cytoscape?

At the very least, you must provide Cytoscape with a list of nodes. You can also provide a list of edges and attributes. You can enter your graph manually or import it from your own files – Cytoscape supports a wide variety of file formats, including text and Excel. Cytoscape can also load pre-configured networks (e.g., from BioGRID) so you can leverage them without having to enter them yourself.

Can Cytoscape Help Me Get Started?

Yes! The Cytoscape team has provided many resources tailored for all kinds of users (see Section 8). To start, there is an Installation Guide and numerous tutorials. Once you load your graph, you may find additional useful analytics and visualizations at the App Store (see Section 6).

Which Cytoscape Version Should I Use?

Generally, you should use v3.1.1 (the newest version). v2.x has been in the field for several years and will be supported into the future, but over time, v3.x and successors will be better supported.

Note that if you have already installed v3.0.0, v3.0.1, v3.0.2, or 3.1, you can install v3.1.1 without uninstalling your previous version.

How to Download and Install v3.x

- Surf to <http://www.cytoscape.org>, click on the Download Cytoscape button, and follow the prompts
- ***Be sure to consult with the Release Notes (viewable via the Release Notes button) for issues and solutions pertaining to your installation***
- Surf to <http://www.cytoscape.org>, click on the Release Notes button, click on a User Manual link, and follow the instructions in the *Launching Cytoscape* section.

2. For Casual Cytoscape Users

If you're confident of your Cytoscape skills, please read the Power Users information (Section 3).

If you're Cytoscape skills are rusty, please continue with the New Users information (Section 1).

3. For Cytoscape Power Users

Is v3.x for Me?

The answer is most likely yes! Most modern v2.x plugins have been converted to v3.x apps, and new v3.x apps are arriving in Cytoscape's App Store far more quickly than v2.x plugins (-- *app* is the v3.x name for the installable features known as a *plugin* in v2.x). To take stock, you can browse the App Store directly (Section 6). Note that some v2.x plugins have been incorporated directly into Cytoscape and need not be downloaded or installed. Consult the App Store to see a list of these apps.

Note that v2.x plugins are not compatible with v3.x apps, and vice versa. If you rely on v2.x plugins that are not yet available for v3.x, you may need to remain with v2.x instead of upgrading to v3.x. Note, also, that the Cytoscape community is converting v2.x plugins to v3.x apps -- you might occasionally revisit the App Store to see if an app you need becomes available.

What's New in v3.x?

As described in the Welcome letter (above), v3.x incorporates a number of new features (Section 7) and productivity enhancing support offerings (Section 8), including the Cytoscape App Store (Section 6). The General Background section (Section 5) describes the rationale, advantages, and issues with v3.x.

Will Support for v2.x Continue?

Cytoscape v2.x will be supported by the Cytoscape community for now, but with fewer upgrades and support options over time.

Can I use v2.x and v3.x at the Same Time?

You can experiment with v3.x by installing it on your workstation – you can use both v2.x and v3.x on the same workstation.

How to Prepare for v3.x

1. If you have already installed a v3.2 beta version, *please uninstall it before proceeding*. The Cytoscape v3.0 and v3.1 versions can coexist with v3.2, and do not need to be uninstalled.
2. Consult the App Store (Section 6) to verify that the apps you need are available. Some apps have become part of the Cytoscape core, and you won't need to download them. Other apps are available in the Cytoscape App store.
3. Browse the New Features table (Section 7) to see what we have added since v2.x.
4. Browse the Future Plans table (Section 9) and Support Offering table (Section 8) to see where the Cytoscape community is going.

How to Download and Install v3.x

- Surf to <http://www.cytoscape.org>, click on the Download Cytoscape button, and follow the prompts
- ***Be sure to consult with the Release Notes (viewable via the Release Notes button) for issues and solutions pertaining to your installation***
- Surf to <http://www.cytoscape.org>, click on the Release Notes button, click on a User Manual link, and follow the instructions in the *Launching Cytoscape* section.

4. For Cytoscape App Developers

Cytoscape manages graph storage and display, and contains a number of common analytic and visualization functions. App Developers may write *apps* (formerly known as *plugins*) to perform problem-specific analytics, graph layouts, or other visualizations.

Before starting app development, you should be familiar with overall Cytoscape usage and functionality – see the Power Users section (Section 3) to start working with v3.x if you have not already done so.

Note that v2.x plugins do not work with v3.x. In order to get the benefits of improved stability, the Cytoscape team found it necessary to change the relationship between plugins/apps and the Cytoscape core. Many plugins have been ported to v3.x (see Section 6). The porting exercise is substantial, and once a plugin is converted to an app, it will continue to work at least until v4.0 is released – there is no v4.0 planned at this time. For more on this topic, see:

http://wiki.cytoscape.org/Cytoscape_3/CoreDevelopment/Motivation

Should I Convert My v2.x Plugin to a v3.x App?

Yes! Cytoscape v2.x will be supported by the Cytoscape community for now, but with fewer upgrades and support options over time.

How Do I Get Started With App Development?

To get started with app development:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper

To convert a v2.x plugin to a v3.x app, refer to:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/PluginPortingGuide

or

<http://f1000research.com/articles/3-150/v1>

For a broader view of the Cytoscape system (including the Core), refer to:

http://wiki.cytoscape.org/Cytoscape_3/CoreDevelopment

Are there App Samples?

Yes! See:

http://opentutorials.cgl.ucsf.edu/index.php/Tutorial:Creating_a_Simple_Cytoscape_3_App

Also note the App Cookbook:

http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/Cytoscape_3_App_Cookbook

5. General Background

The overall mission of Cytoscape is to be a freely available worldwide asset supporting network analysis and visualization for systems biology science.

The major focus of v3.x is the modularization and rationalization of code to solve stability issues in v2.x encountered as multiple developers pursued multiple agendas. Under v2.x, internal programmatic interfaces evolved from one release to the next, leading to the failure of working plugins over time and negative interactions between otherwise working plugins. Ultimately, this resulted in loss of programmer and user productivity, and undermined community confidence in Cytoscape.

v3.x addresses these issues by adopting modular coding practices promoted by the OSGi architectural framework¹. This enables both the Cytoscape core and externally developed apps (formerly called *plugins*) to evolve independently without compromising unrelated functionality. At the logical level, Cytoscape leverages OSGi precepts to produce v3.x APIs having cleaner and clearer demarcations between functional areas. At the deployment level, OSGi enables on-the-fly substitution of one processing element for another (e.g., apps) in order to tailor Cytoscape to meet user requirements at runtime without reinstalling or reconfiguring Cytoscape.

Creating v3.x occupied an international team of engineers for over a year, and represents a strong investment toward reducing future development and support costs, and increasing reliability and evolvability. With the completion of v3.x, we expect to leverage v3.x as a platform to satisfy the evolving needs of multiple stakeholder groups, and as a platform enabling research on leading edge analysis and visualization techniques. v3.x is the intended successor to v2.x, with development and support of v2.x expected to diminish and disappear over time in favor of v3.x.

We believe that users will benefit most directly from the v3.x by:

- experiencing fewer core and app bugs from one release to the next
- the availability of more and richer apps (due to developers spending less time tracking and fixing bugs)
- more core features with higher biological and logistical value (due to improved flexibility provided by interface-driven development)

During the creation of v3.x, a number of features were added, as listed in the New Feature Table below.

¹ www.osgi.org – also used as the basic framework for Eclipse and numerous commercial products

6. The Cytoscape App Store

The Cytoscape App Store is new in v3.x and targets broad demographics through specific features:

Cytoscape User	Benefit
New or Casual	Find Cytoscape functionality easier via familiar paradigms of search and tag navigation
Power or Frequent	Find and install apps easier and provide crowdsourced rating and reviews
App Developers	Have a dedicated presence for their work, including a direct line to their users and critical statistics to justify their development on the Cytoscape platform.

Furthermore, the App Store supports contests and social incentives that encourage critical technical and community development, including porting, analytics, rating/reviewing, and so on.

Note that the App Store supports only Cytoscape v3.x and does not contain plugins compatible with earlier Cytoscape versions. Likewise, plugins that work with earlier versions are not compatible with v3.x – the App Store (<http://apps.cytoscape.org>) contains a list of apps available for v3.x.

To access the App Store from Cytoscape, access Cytoscape's new Apps menu and choose the App Manager menu item.

7. New Feature Table

Feature	Beneficiary	Manual Section
R/Matlab/Python REST integration	Power users	http://apps.cytoscape.org/apps/cyrest
Fast network loads	All users	
Publish layouts to web	All users	23. Publish Your Data
Create entire web sites	Power users	23. Publish Your Data
Fast network merge	All users	
Java 8 support (experimental)	All users	3. Launching Cytoscape
Fast startup times	All users	
Nodes displayed as charts	Power users	14.5.6. Creating Node Charts
Load larger networks	Power users	3. Launching Cytoscape

8. Support Offering Table

Feature	Purpose	Beneficiary	URL
Cytoscape Google Groups Helpdesk Discuss	Community response on best effort basis	All users	mailto:cytoscape-helpdesk@googlegroups.com mailto:cytoscape-discuss@googlegroups.com
User Manual	Comprehensive manual	All users	http://wiki.cytoscape.org/Cytoscape_3/UserManual
Installation Guide	Cytoscape installation	New/Casual users	Section 3 of http://wiki.cytoscape.org/Cytoscape_3/UserManual
Tutorials	Streamlined productivity	All users	http://opentutorials.cgl.ucsf.edu/index.php/Portal:Cytoscape3
Presentation Slides on Web Site	Enable casual learning	All users	https://speakerdeck.com/cytoscape
Code Snippets	Commonly used App development code	App Developers	http://wiki.cytoscape.org/Cytoscape_3/AppDeveloper/Cytoscape_3_App_Cookbook

9. Future Plans Table

Feature	Purpose	Beneficiary	Timeframe
R/Mathlab/Python integration	Enable model access, analysis, and visualization from popular tools	Power users	Spring 2015 – Python available NOW
Large network handling	Improved layout, interactive, and load performance	Power users	Spring 2015 – partially available now
Movies	Clear and interactive engagement	New/Casual users	TBD
Additional Visualizations			TBD
Additional Analytics			TBD
Wallpaper Collection	Publicize good results associated with Cytoscape	All users, potential users, general public	TBD
Case Studies	Streamlined productivity	New/Casual Users	TBD
Cytoscape Clinics New/Casual User Power User App Developer	Chat rooms give instant response based on community availability	All users	TBD

10. Issues Addressed

The following issues were addressed in v3.2. Details for each issue can be found at Cytoscape's RedMine issue repository:

<http://code.cytoscape.org/redmine>

Issue	Description
1648	Improve interactive reponse for large networks (and other cool features...)
1788	Race condition sometimes causes apps to initialize before core bundles
1887	Remove unused image files in swing-application-impl
2143	Need to support unattended install
2192	Increase Startup Speed of Cytoscape 3
2196	Destroying single view after creating multiple CyNetworkView causes collection title to be red
2243	Table Browser does not display third-party network tables
2258	Create a Grunt Project to display exported Cytoscape.js JSON files
2373	No way to permanently import style
2374	Network View's style not updated after deleting a mapped column
2410	Table export/merge greyed out until table selector clicked
2415	Need a function to export session as HTML5
2426	Add "required" flag to Tunable annotations
2434	PSI-MI file writer creates invalid PSIMI 2.5 XML file
2440	Execution Icon in Launcher is Java rather than Cytoscape in Ubuntu
2493	Merge networks needs to provide better progress
2499	Cut option doesn't immediately remove node/edge
2505	Style options should be ordered
2517	Creating an OSGi Bundle tutorial needs to be updated.
2524	Add support for importing attribute files without .attrs extension
2525	Problems with annotations impl
2530	Rendering Engine ID of network views not saved in the session file
2532	Install From File shows duplicate app entries in App Manager Currently Installed
2538	Merging large networks needs better user feedback
2572	Cytoscape fails to apply edge style when loading session
2583	Command line dialog does not display status messages during task execution
2584	Importing table or network from command line does not work if delimiter is COMMA
2594	Issue using group through command line dialog - collapse command doesn't show ...
2597	Commands do not check input parameters before running tasks
2600	[Table import] Changing the list delimiter for a column causes the primary key mapping ...
2626	Add Node Charts to VizMapper
2644	Exceptions thrown by ContextMenuFactories cause node/edge/network context menus ...
2647	Once a group is created and collapsed loose visual style on the edges
2671	Edge Color is not refreshing
2685	Data-type error while importing tables from file
2701	CySessionManagerImpl.restoreProperties() re-registers CyProperty services
2721	-S/--script command line option does not handle file paths correctly
2726	Broken network view when creating subnetwork ("__Annotations column already exists" error)
2727	Java 8 support

2738 IllegalArgumentException when running Edit>Remove Duplicate Edges with column create ...
2749 command line export view no longer working
2761 Image annotations disappear on zoom
2762 Shape annotations (including images) should support a "Shift" modifier
2763 Network Columns of type List of Longs with the name ending in ".SUID" don't get SUIDs ...
2764 "help table export" command
2767 Table Import from File half-dies when parsing fails
2770 Node positions aren't getting updated as part of an undo
2775 Format string argument for tunable numbers
2777 Font face dialog not appearing
2781 pre-loaded edgelinkouturls don't work
2785 Annotation setFontSize broken
2786 Relying on floating-point numbers causing inconsistent rounding behavior
2788 Edge transparency setting in style has no effect
2793 Annotation setFillOpacity has no effect
2799 CyNetworkView updateView not finding nodes
2802 JSON Exporter does not export color passthrough mapping
2804 Command Line Error Message
2805 Import GO Annotation feature does not work in 3.1.1
2807 Missing style buttons
2808 JSON Importer needs smarter number value parser
2811 ArrayIndexOutOfBoundsException when rendering INVISIBLE elements
2824 Cannot add networks to existing Network Collection when loading JSON
2835 Collapsing a group that contains another expanded group does not work sometimes
2836 Update embedded resource file to Cytoscape.js 2.3.x
2839 Child group node of another collapsed group incorrectly displayed after loading the session
2841 Ding throws concurrent modification exception
2842 Ding throws null pointer exception
2843 Cytoscape falls into deadlock when selecting nodes/edges right after zoom/pan
2844 Can't uninstall an app that was just installed on Windows without restarting Cytoscape
2847 Welcome screen on Linux has inconsistent color scheme...
2848 Deleted nodes/edges that were selected always displayed in table browser
2849 Collapsing some groups causes an exception to be thrown on 3.2.0-SNAPSHOT
2851 Collapsing group is 10k times slower than 3.1.x
2852 Bypass values set to group elements lost after saving and reopening a session
2853 Collapsing/Expanding Groups throws series of exceptions
2854 Can't scroll the command line dialog up
2855 setNetworkPointer error (NODE_NESTED_NETWORK_IMAGE_VISIBLE VisualProperty)
2856 Dead link on website
2857 Ding throws NPE at random
2858 Chart implementation backtraces when zooming past high level of detail
2860 Need a warning message to users when they try to export large networks into HTML5 session
2863 Default Custom Graphics does not work in Cytoscape 3.2
2864 cytoscape:script in OSGi Console is not available any longer
2866 Cannot save session with visual style containing charts

2868 Create better sample session file using latest features
2869 Layout settings are no longer working
2876 Custom Graphics Passthrough Mapping crashes Cytoscape
2879 Need 32-bit support for 3.2
2881 Edge Weighted layout problems
2886 Export to Simple HTML5 exports first network, not the selected one
2887 Continuous mapping exporter for JSON Style does not always work
2888 Gradient layer does not work
2889 No way to perform weighted layouts using command line.
2891 Issues saving with multiple networks but error only comes up sometimes
2899 Images are not saved in session file in some situations